



Autoclave Biological Indicator Test Log

Autoclave Identification:				Autoclave Location:		Name Of PI/Dept:	
S/N of Autoclave				Building Name		Dr. Murray	
Date of Test	Name of Cycle	Duration Of Cycle	Time	Temp.	Pressure	Result of Test +/- (Color)	Initials
8/11/25	Waste Treat Cycle	45 min.	10 a.m.	121° C (250° F)	30 psi	+ Yellow	JM
9/15/25	Waste Treat Cycle	45 min.	10 a.m.	121° C (250° F)	30 psi	- Purple	JM
10/13/25	Waste Treat Cycle	45 min.	10 a.m.	121° C (250° F)	30 psi	- Purple	JM

EXAMPLE

