



Biological Safety Cabinet Use, Maintenance, and Disinfection

Building/Room Number: _____ BSC Make/Model: _____ MEF#: _____

Approved Disinfectants: 1. Clorox Bleach Germicidal Cleaner 2. CaviCide 3. Sporicidal 4. Steri-Perox 5. Decon-Ahol (70% IPA)

Notes:

- Cleaning by Clorox bleach Germicidal Cleaner needs to be followed by a Decon-Ahol cleaning after the contact time has been met to remove residue.
- Rotate the use of at least two approved disinfectants sing one per day of use of equipment.

Date	Time In	Time Out	Disinfectant Used	Product Used	Performed By	Activity/ Comments
9/18/25	9:00am	10:00am	CaviCide	Ethanol (70% IPA)	Dr. Murray	Detailed Cleaning
9/25/25	9:00am	10:00am	CaviCide	Ethanol (70% IPA)	Dr. Murray	Detailed Cleaning

EXAMPLE

Date	Time In	Time Out	Disinfectant Used	Product Used	Performed By	Activity/ Comments



Environmental Health & Safety