



Environmental Health & Safety

Biological Centrifuge Use, Maintenance, and Disinfection

Building/Room Number: _____ Make/Model: _____ MEF#: _____

Approved Disinfectants: 1. Clorox Bleach Germicidal Cleaner 2. CaviCide 3. Decon-Ahol (70% IPA)

Notes:

- Be sure to always inspect the rotors and the rotor chambers for damage.
- Lubrication of the buckets, pivots and rubber seals are crucial to maintain an operational centrifuge.

Date	Time In	Time Out	Disinfectant Used	Performed By	Activity/ Comments
9/18/25	9:00am	10:00am	CaviCide; Ethanol (70% IPA)	Dr. Murray	Cleaned & Oiled
9/25/25	9:00am	10:00am	CaviCide; Ethanol (70% IPA)	Dr. Murray	Cleaned & Oiled

Centrifuge Cleaning Log
Revision Date: 10/6/2025

www.riskmanagement.edu
940-565-2109

Please insert this sheet of paper into a transparent document sleeve.

Version 1.0.0

Date	Time In	Time Out	Disinfectant Used	Performed By	Activity/ Comments