



Environmental Health & Safety

Biological Incubator Use, Maintenance, and Disinfection

Building/Room Number: _____ Make/Model: _____ MEF#: _____

Approved Disinfectants: 1. CaviCide 2. Decon-Ahol (70% IPA)

Notes:

- Cleaning with soapy water needs to be followed by a Quaternary ammonium disinfectant.
- After the Quaternary ammonium contact time has been met, the residue must be removed with 70% ethanol

Date	Time In	Time Out	Disinfectant Used	Performed By	Activity/ Comments
9/18/25	9:00am	10:00am	CaviCide; Ethanol (70% IPA)	Dr. Murray	Cleaned & Filled w/ Water
9/25/25	9:00am	10:00am	CaviCide; Ethanol (70% IPA)	Dr. Murray	Cleaned & Filled w/ Water

EXAMPLE

