

UNIVERSITY OF NORTH TEXAS AUTOCLAVE USE LOG

Autoclave ID#: _____ Building: _____ Room#: _____

• **Weight:** The total weight of the load is considered waste if waste is any part of a load.

• **Minimum Parameters for Solid Waste:** Pressure-->15 psig Temperature-->123°C Time-->70 minutes autoclave time

| Date of Treatment | Principal Investigator (PI) | Name of Person Treating Waste (Printed) | Initials | Check all that Apply: | NO WASTE | Plant waste ONLY | Microbio Waste | Sharps Waste | Pathological Waste | Load Sterility verified by CI ? | Weight of waste load | LIQ | SOL | Pre-vac | Run Time (min)/ Temp (°C) |
|-------------------|-----------------------------|---|----------|-----------------------|----------|------------------|----------------|--------------|--------------------|---------------------------------|----------------------|-----|-----|---------|---------------------------|
| SAMPLE | Shoff | M. Shoff | MS | | | | ✓ | | | Yes | 20 lbs | | ✓ | | 70 min/123C |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |

NOTE: Autoclave Waste Treatment Testing Records are to be kept three (3) years for review by the Texas Commission on Environmental Quality at any time.

Forward this completed form to the BSO at biosafety@unt.edu.